


Product Name: Recombinant Human FCRL4 Protein
Catalog #: PEH5989



Summary

Name	FCRL4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96PJ5
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LGDETRLPPA PPGGESSHSI CPAQVELQSL YVDVHPKKGD LVYSEIQTQ LGEEEEANTS RTLLEDKDVS VVYSEVKTQH PDNSAGKISS

image



Background

Alternative Names Fc receptor-like protein 4 (FcR-like protein 4) (FCRL4) (Fc receptor homolog 4) (FcRH4) (IFGP family protein 2) (hIFGP2) (Immune receptor translocation-associated protein 1) (CD antigen CD307d)

Note

For Research Use Only , Not for Diagnostic Use