


Product Name: Recombinant Human MED24 Protein
Catalog #: PEH5998



Summary

Name	MED24_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O75448
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ADQRCNCDCT NFLLQECGKQ GLLSEASVNN LMAKRKADRE HAPQQKSGEN ANIQPNLIQLI LRAEPTVTNI LKTMDADHSK

image



Background

Alternative Names	Mediator of RNA polymerase II transcription subunit 24 (Activator-recruited cofactor 100 kDa component) (ARC100) (Cofactor required for Sp1 transcriptional activation subunit 4) (CRSP complex subunit 4) (Mediator complex subunit 24) (Thyroid hormone receptor-associated protein 4) (Thyroid hormone receptor-associated protein complex 100 kDa component) (Trap100) (hTRAP100) (Vitamin D3 receptor-interacting protein complex 100 kDa component) (DRIP100)
--------------------------	---

Product Name: Recombinant Human MED24 Protein
Catalog #: PEH5998



Note

For Research Use Only , Not for Diagnostic Use