


Product Name: Recombinant Human CA4 Protein
Catalog #: PEH5999



Summary

Name	CAH4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & N-His
Accession #	P22748
Host	E.coli
Species	Human
Molecular Mass	42 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SHWCYEVQAE SSNYPCLVPV KWGGNCQKDR QSPINIVTTK AKVDKKLGRF FFSGYDKKQT WTVQNNGHSV MMLLENKASI SGGGLPAPYQ AKQLHLHWSD LPYKGSEHSL DGEHFAMEMH IVHEKEKGT S RNVKEAQDPE DEIAVLAFLV EAGTQVNEGF QPLVEALSNI PKPEMSTTMA ESSLLDLLPK EEKLRHYFRY LGSLTTPTCD EKVVWTVFRE PIQLHREQIL AFSQKLYYDK EQTVSMKDNV RPLQQLGQRT

image



Background

Alternative Names	Carbonic anhydrase 4 (EC 4.2.1.1) (Carbonate dehydratase IV) (Carbonic anhydrase IV) (CA-IV)
--------------------------	--

Product Name: Recombinant Human CA4 Protein
Catalog #: PEH5999



Note

For Research Use Only , Not for Diagnostic Use