


Product Name: Recombinant Human LGALS14 Protein
Catalog #: PEH6005



Summary

Name	PPL13_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8TCE9
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MSSLPVPYTL PVSLPVGSCV IITGTPILTF VKDPQLEVNF YTGMDSDI AFQFRLHFGH PAIMNSCVFG IWRYEEKCY Y LPFEDGKPF E LCIYVRHKEY KVMVNGQRIY NFAHRFP PAS VKMLQVFRDI SLTRVLISD

image



Background

Alternative Names Placental protein 13-like (Charcot-Leyden crystal protein 2) (CLC2) (Galectin-14) (Gal-14)

Note

For Research Use Only , Not for Diagnostic Use