


Product Name: Recombinant Human P3H1 Protein
Catalog #: PEH6009



Summary

Name	P3H1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q32P28
Host	E.coli
Species	Human
Molecular Mass	14 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DSWTPEEVIP KRLQEKQKSE RETAVRISQE IGNLMKEIET LVEEKTRESL DVSRLTREGG PLLYEGISLT MNSKLLNGSQ RVVMDGVISD HECQELQRLT NVAATSGDGY

image



Background

Alternative Names	Prolyl 3-hydroxylase 1 (EC 1.14.11.7) (Growth suppressor 1) (Leucine- and proline-enriched proteoglycan 1) (Leprecan-1)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use