


Product Name: Recombinant Human MED12L Protein
Catalog #: PEH6011



Summary

Name	MD12L_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q86YW9
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IKYEEQHLL LYHTHPMPKP RSYLQPLPL PEEEEEEPT SPVSQEPERK SAELSDQGKT TTDEEKTKG

image



Background

Alternative Names	Mediator of RNA polymerase II transcription subunit 12-like protein (Mediator complex subunit 12-like protein) (Thyroid hormone receptor-associated-like protein) (Trinucleotide repeat-containing gene 11 protein-like)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human MED12L Protein
Catalog #: PEH6011

