

Product Name: Recombinant Human POU3F3 Protein
Catalog #: PEH6015



Summary

Name	PO3F3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P20264
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	HPPHPHHAQG PPHHGGGGGG AGPGLNSHDP HSEDTPTSD DLEQFAKQFK QRRIKLGFTQ ADVGLALGTL YGNVFSQTTI CRFEALQLSF KNMCKLKPLL NKWLEEADSS TGSPTSIDKI AAQGRKRKRR TSIEVSVKGA LESHFLKCPK PSAQEITNLA DSLQLEKEVV RVWFCNRRQK EKRMTPPGIQ QQTPDDVYSQ VGTVSADTPP

image

SDS

Background

Alternative Names	POU domain, class 3, transcription factor 3 (Brain-specific homeobox/POU domain protein 1) (Brain-1) (Brn-1) (Octamer-binding protein 8) (Oct-8) (Octamer-binding transcription factor 8) (OTF-8)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use