


Product Name: Recombinant Human LPCAT1 Protein
Catalog #: PEH6019



Summary

Name	PCAT1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8NF37
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DLDRYSERAR MKGGEEKIGIA EFAASLEVPV SDLLEDMFSL FDESGSGEVD LRECVVALSV VCRPARTLDT IQLAFKMYGA QEDGSVGEGD

image



Background

Alternative Names	Lysophosphatidylcholine acyltransferase 1 (LPC acyltransferase 1) (LPCAT-1) (LysoPC acyltransferase 1) (EC 2.3.1.23) (1-acylglycerol-3-phosphate O-acyltransferase) (EC 2.3.1.51) (1-acylglycerophosphocholine O-acyltransferase) (1-alkenylglycerophosphocholine O-acyltransferase) (EC 2.3.1.25) (1-alkylglycerophosphocholine O-acetyltransferase) (EC 2.3.1.67) (Acetyl-CoA:lyso-platelet-activating factor acetyltransferase) (Acetyl-CoA:lyso-PAF acetyltransferase) (Lyso-PAF acetyltransferase) (LysoPAFAT)
--------------------------	---

Product Name: Recombinant Human LPCAT1 Protein
Catalog #: PEH6019



(Acyltransferase-like 2) (Phosphonoformate immuno-associated protein 3)

Note

For Research Use Only , Not for Diagnostic Use