


Product Name: Recombinant Human PF4V1 Protein
Catalog #: PEH6021



Summary

Name	PF4V_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P10720
Host	E.coli
Species	Human
Molecular Mass	9 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FARAEAEEDG DLQCLCVKTT SQVRPRHITS LEVIKAGPHC PTAQLIATLK NGRKICLDLQ ALLYKKIIKE HLES

image



Background

Alternative Names	Platelet factor 4 variant (C-X-C motif chemokine 4 variant) (CXCL4L1) (PF4alt) (PF4var1) [Cleaved into: Platelet factor 4 variant(4-74); Platelet factor 4 variant(5-74); Platelet factor 4 variant(6-74)]
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use