



Product Name: Recombinant Human LGALS9 Protein
Catalog #: PEH6027



Summary

Name	LEG9_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	O00182
Host	E.coli
Species	Human
Molecular Mass	52 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MAFSGSQAPY LSPAVPFSGT IQGGLQDGLQ ITVNGTVLSS SGTRFAVNFQ TGFSGNDIAF HFNPRFEDGG YVVCNTRQNG SWGPEERKTH MPFQKGMPFD LCFLVQSSDF KVMVNGILFV QYFHRVPFHR VDTISVNGSV QLSYISFQNP RTVPVQPAFS TVPFSQPVCF PPRPRGRRQK PPGVWPANPA PITQTVIHTV QSAPGQMFST PAIPPMMPYPH PAYPMPFITT ILGGLYPSKS ILLSGTVLPS AQRFHINLCS GNHIAFHLNP RFDENAVVRN TQIDNSWGSE ERS�PRKMPF VRGQSFSVWI LCEAHCLKVA VDGQHLFEYY HRLRNLPTIN RLEVGGDIQL THVQT

image



Background

Alternative Names	Galectin-9 (Gal-9) (Ecalectin) (Tumor antigen HOM-HD-21)
--------------------------	--

Product Name: Recombinant Human LGALS9 Protein
Catalog #: PEH6027



Note

For Research Use Only , Not for Diagnostic Use