


Product Name: Recombinant Human METTL5 Protein
Catalog #: PEH6028



Summary

Name	METL5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NRN9
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ENKVVADLGC GCGVLSIGTA MLGAGLCVGF DIDEALEIF NRNAEEFELT NIDMVQCDVC LLSNRMSKSF DTVIMNPPFG TKNNKGTDMA FLKTALEMAR TAVYSLHKSS TREHVQKKAA EWKIKIDIIA ELRYDLPASY KFHKKKSVDI

image



Background

Alternative Names	rRNA N6-adenosine-methyltransferase (Methyltransferase-like protein 5)	METTL5	(EC 2.1.1.-)
--------------------------	---	--------	--------------

Note

For Research Use Only , Not for Diagnostic Use