



Product Name: Recombinant Human POLR2J Protein
Catalog #: PEH6031



Summary

Name	RPB11_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	P52435
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MNAPPAFESF LLFEGEKKIT INKDTKVPNA CLFTINKEDH TLGNIKSQL LKDPQVLFAG YKVPHPLEHK IIRVQTTPD YSPQEAFTNA ITDLISELSL LEERFRVAIK DKQEGIE

image



Background

Alternative Names	DNA-directed RNA polymerase II subunit RPB11-a (RNA polymerase II subunit B11-a) (RPB11a) (DNA-directed RNA polymerase II subunit J-1) (RNA polymerase II 13.3 kDa subunit)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human POLR2J Protein
Catalog #: PEH6031

