

Product Name: Recombinant Human RGS16 Protein
Catalog #: PEH6032



Summary

Name	RGS16_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O15492
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TRLGIFLHKS ELGCDTGSTG KFEWGSKHSK ENRNFSEDVL GWRESFDLLL SSKNGVAAFH AFLKTEFSEE NLEFWLACEE FKKIRSATKL ASRAHQIFEE FICSEAPKEV NIDHETHELT RMNLQTATAT CFDAAQGKTR TLMEKDSYPR FLKSPAYRDL AAQASAASAT LSSCSLDEPS

image



Background

Alternative Names Regulator of G-protein signaling 16 (RGS16) (A28-RGS14P) (Retinal-specific RGS) (RGS-r) (hRGS-r) (Retinally abundant regulator of G-protein signaling)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human RGS16 Protein
Catalog #: PEH6032

