


Product Name: Recombinant Human PRTN3 Protein
Catalog #: PEH6040



Summary

Name	PRTN3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	P24158
Host	E.coli
Species	Human
Molecular Mass	37 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AEIVGGHEAQ PHSRPYMASL QMRGNPGSHF CGGTLIHPSF VLTAAHCLRD IPQRLVNVVL GAHNVRTQEP TQQHFSVAQV FLNNYDAENK LNDVLLIQLS SPANLSASVA TVQLPQQDQP VPHGTQCLAM GWGRVGAHDP PAQVLQELNV TVVTFPCRPH NICTFVPRRK AGICFGDSGG PLICDGIIQG IDSFVIWGCA TRLFPDFTR VALYVDWIRS TLRR

image



Background

Alternative Names Myeloblastin (EC 3.4.21.76) (AGP7) (C-ANCA antigen) (Leukocyte proteinase 3) (PR-3) (PR3) (Neutrophil proteinase 4) (NP-4) (P29) (Wegener autoantigen)

Note

Product Name: Recombinant Human PRTN3 Protein
Catalog #: PEH6040



For Research Use Only , Not for Diagnostic Use