


Product Name: Recombinant Human LTK Protein
Catalog #: PEH6045



Summary

Name	LTK_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	P29376
Host	E.coli
Species	Human
Molecular Mass	31 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AGTSVVVTVG AAGQLRGVQL WRVPGPGQYL ISAYGAAGGK GAKNHLSRAH GVFVSAIFSL GLGESLYILV GQQGEDACPG GSPESQLVCL GESRAVEEHA AMDGSEGVPV SRRWAGGGGG GGGATYVFRV RAGELEPLLV AAGGGGRAYL RPRDRGRTQA SPEKLENRSE

image



Background

Alternative Names Leukocyte tyrosine kinase receptor (EC 2.7.10.1) (Protein tyrosine kinase 1)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human LTK Protein
Catalog #: PEH6045

