

Product Name: Recombinant Human NKX1-2 Protein
Catalog #: PEH6046



Summary

Name	NKX12_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UD57
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SRDPVRQLET PDAAGPGAGQ ASPLEGSEAE EEEDAEDPRR PRLRERAARL LPGLARSPDA PAGALASGEP CEDGGGGPVR SPPGSPGSPR PRRRRLEPNC AKPRRARTAF TYEQLVALEN KFRATRYLSV CERLNLALSL SLTETQVKIW FQNRRTKWKK QNPGADGAAQ VGGGAPQPGA AGGGGGGGSG

image

SDS

Background

Alternative Names	NK1 transcription factor-related protein 2 (Homeobox protein SAX-1) (NKX-1.1)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use