


Product Name: Recombinant Human LDB1 Protein
Catalog #: PEH6050



Summary

Name	LDB1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q86U70
Host	E.coli
Species	Human
Molecular Mass	14 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KEPPNGNAFP PFHPGTMLDR DVGPTPMYPP TYLEPGIGRH TPYGNQTDYR IFELNKRLQN WTEECDNLWW DAFTEFFED DAMLTITFCL EDGPKRYTIG

image



Background

Alternative Names LIM domain-binding protein 1 (LDB-1) (Carboxyl-terminal LIM domain-binding protein 2) (CLIM-2) (LIM domain-binding factor CLIM2) (hLdb1) (Nuclear LIM interactor)

Note

For Research Use Only , Not for Diagnostic Use