


Summary

Name	LIRA5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	A6NI73
Host	E.coli
Species	Human
Molecular Mass	38 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AEPGSVISRG NSVTIRCQGT LEAQEYRLVK EGSPEPWDTQ NPLEPKNKAR FSIPSMTEHH AGRYRCYYYS PAGWSEPSDP LELVVTGFYN KPTLSALPSP VVTSGENVTL QCGSRLRFDR FILTEEGDHK LSWTLDSQLT PSGQFQALFP VGPVTPSHRW MLRCYGSRRH ILQVWSEPSD LLEIPVSGAA DNLSPSQNKs DSGTASHLQD

image



Background

Alternative Names	Leukocyte immunoglobulin-like receptor subfamily A member 5 (CD85 antigen-like family member F) (Immunoglobulin-like transcript 11) (ILT-11) (Leukocyte immunoglobulin-like receptor 9) (LIR-9) (CD antigen CD85f)
--------------------------	--

Product Name: Recombinant Human LILRA5 Protein
Catalog #: PEH6052



Note

For Research Use Only , Not for Diagnostic Use