


Product Name: Recombinant Human POLR3E Protein
Catalog #: PEH6054



Summary

Name	RPC5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NVU0
Host	E.coli
Species	Human
Molecular Mass	14 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KHPDVVQRQH MLWTGIQAKL EKVYNLVKET MPKKPDAQSG PAGLVCGDQR IQVAKTKAQQ NHALLERELQ RRKEQLRVPA VPPGVRIKEE PVSEEGEED EQEAEEPPMD

image



Background

Alternative Names DNA-directed RNA polymerase III subunit RPC5 (RNA polymerase III subunit C5) (DNA-directed RNA polymerase III 80 kDa polypeptide)

Note

For Research Use Only , Not for Diagnostic Use