


Product Name: Recombinant Human POLDIP2 Protein
Catalog #: PEH6056



Summary

Name	PDIP2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y2S7
Host	E.coli
Species	Human
Molecular Mass	32 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RLYDRDVASA APEKAENPAG HGSKEVKGKT HTYYQVLIDA RDCPHISQRS QTEAVTFLAN HDDSRALYAI PGLDYVSHED ILPYTSTDQV PIQHELPERF LLYDQTKAPP FVARETLRAW QEKNHPWLEL SDVHRETTEN IRVTVIPFYM GMREAQNSHV YWWRYCIRLE NLDSVVQLR ERHWRFISLS GTLETVRGRG VVGREPVLSK EQPAFYSSH VSLQASSGHM WGTFRFERPD GSHFDVRIPP FSLESNKDEK

image



Background

Alternative Names Polymerase delta-interacting protein 2 (38 kDa DNA polymerase delta interaction protein) (p38)

Product Name: Recombinant Human POLDIP2 Protein
Catalog #: PEH6056



Note

For Research Use Only , Not for Diagnostic Use