


**Product Name: Recombinant Human CYC1 Protein**  
**Catalog #: PEH6058**



## Summary

<b>Name</b>	CY1_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	N-SUMO & C-His
<b>Accession #</b>	P08574
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	26 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	PPSYPWSHRG LLSSLDHTSI RRGFQVYKQV CASCHSMDFV AYRHLVGVCY TEDEAKELAA EVEVQDGPNE DGEMFMRPGK LFDYFPKPYP NSEARAANN GALPPDLSYI

## image



## Background

<b>Alternative Names</b>	Cytochrome c1, heme protein, mitochondrial (EC 7.1.1.8) (Complex III subunit 4) (Complex III subunit IV) (Cytochrome b-c1 complex subunit 4) (Ubiquinol-cytochrome-c reductase complex cytochrome c1 subunit) (Cytochrome c-1)
--------------------------	--

## Note

**Product Name: Recombinant Human CYC1 Protein**  
**Catalog #: PEH6058**



---

For Research Use Only , Not for Diagnostic Use