


Product Name: Recombinant Human TMEM161B Protein
Catalog #: PEH6063



Summary

Name	T161B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q8NDZ6
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QHPTEELRI LAGKQQKGKT KKDRKYNGHI ESKPLTIPKD IDLHLETKSV TEVDTLALHY

image



Background

Alternative Names Transmembrane protein 161B

Note

For Research Use Only , Not for Diagnostic Use