


Product Name: Recombinant Human SGMS1 Protein
Catalog #: PEH6071



Summary

Name	SMS1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q86VZ5
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KVADWLLENA MPEYCEPLEH FTGQDLINLT QEDFKKPLC RVSSDNGQRL LDMIETLKME HHLEAHKNGH ANGHLNIGVD IPTPDGSFSI KIKPNGMPNG YRKEMIKIPM PELERSQYPM

image



Background

Alternative Names	Phosphatidylcholine:ceramide cholinephosphotransferase 1 (EC 2.7.8.27) (Medulla oblongata-derived protein) (Protein Mob) (Sphingomyelin synthase 1) (Transmembrane protein 23)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human SGMS1 Protein
Catalog #: PEH6071

