


Product Name: Recombinant Human SST Protein
Catalog #: PEH6075



Summary

Name	SMS_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-His
Accession #	P61278
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LRQFLQKSLA AAAGKQELAK YFLAELLSEP NQTENDALEP EDLSQAAEQD EMRLELQRSA NSNPAMAPRE RKAGCKNFFW

image



Background

Alternative Names Somatostatin (Growth hormone release-inhibiting factor) [Cleaved into: Somatostatin-28; Somatostatin-14 (SST-14); Neuronostatin (NST)]

Note

For Research Use Only , Not for Diagnostic Use