


Product Name: Recombinant Human SLC9C1 Protein
Catalog #: PEH6078



Summary

Name	SL9C1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q4G0N8
Host	E.coli
Species	Human
Molecular Mass	33 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LGYLEYDHPE IAVTVKTKEE INVMLNMATE ILKAFGLKGI ISKTEGAGIN KLIMAKKKEV LDSQSIIRPL TVEEVLYHIP WLDKNKDYIN FIQEKAKVVT FDCGNDIFEE GDEPKGIIYI ISGMVKLEKS KPGLGIDQMV ESKEKDFPII DTDYMLSGEI IGEINCLTNE PMKYSATCKT VVETCFIPKT HLYDAFEQCS PLIKQKMWLK LGLAITARKI REHLSYEDWN YNMQKLKLSNI YVVDIPMSTK TDIYDENLIY VILIHGAVED

image



Background

Alternative Names	Sodium/hydrogen exchanger 10 (Na ⁺)/H ⁺ exchanger 10) (NHE-10) (Solute carrier family 9 member 10) (Solute carrier family 9 member C1) (Sperm-specific Na ⁺ /H ⁺ exchanger) (sNHE)
--------------------------	---

Product Name: Recombinant Human SLC9C1 Protein
Catalog #: PEH6078



Note

For Research Use Only , Not for Diagnostic Use