


Product Name: Recombinant Human SLAMF8 Protein
Catalog #: PEH6084



Summary

Name	SLAF8_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q9P0V8
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QPWTQTLQLK VYDAVPRPVV QVFIAVERDA QPSKTCQVFL SCWAPNISEI TYSWRRETTM DFGMEPHSLF TDGQVLSISL GPGDRDVAYS CIVSNPVSWD LATVTPWDSC HHEAAPGKAS

image



Background

Alternative Names	SLAM family member 8 (B-lymphocyte activator macrophage expressed) (BCM-like membrane protein) (CD antigen CD353)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use