


Product Name: Recombinant Human NTMT2 Protein
Catalog #: PEH6090



Summary

Name	NTM1B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q5VVY1
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IELSSPDIQA SQKFLRKFGV GPGRAGTDCA LDCGSGIGRV SKHVLLPVFN SVELVDMMES FLEAQNYLQ VKGDKVESYH CYSLQEFTPP FRRYDVIWIQ WVSGHLTDKD LLAFSLRCRD GLKENGIIIL KDNVAREGCI LDLSDSSVTR DMDILRSLIR KSGLVVLGQE

image



Background

Alternative Names Alpha N-terminal protein methyltransferase 1B (EC 2.1.1.299) (Methyltransferase-like protein 11B) (X-Pro-Lys N-terminal protein methyltransferase 1B) (NTM1B)

Note

Product Name: Recombinant Human NTMT2 Protein
Catalog #: PEH6090



For Research Use Only , Not for Diagnostic Use