



Summary

| | |
|------------------------------|---|
| Name | P0DUX1 · SPD12_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | N-His |
| Accession # | P0DUX1 |
| Host | E.coli |
| Species | Human |
| Molecular Mass | 16 kDa |
| Buffer | Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | SHQPQPQEEQ SPQRSTSGYP LQEVVDDEVS GPSAPGVDPS PPRRSLGWKR KRECLDESDD EPEKELAPEP EETWVAETLC GLKMKAKRRR VSLVLPEYYE AFNRLLEDPV IKRLLAWDKD |

image



Background

Alternative Names Speedy protein E12

Note

For Research Use Only , Not for Diagnostic Use