


Product Name: Recombinant Human SPG7 Protein
Catalog #: PEH6106



Summary

Name	SPG7_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UQ90
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	HEARLLVAKA YRHTEKVLQD NLDKLQALAN ALLEKEVINY EDIEALIGPP PHGPKKMIAP QRWIDAQREK QDLGEEETEE

image



Background

Alternative Names	Paraplegin (EC 3.4.24.-) (Cell matrix adhesion regulator) (Spastic paraplegia 7 protein)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use