


Product Name: Recombinant Human OIP5 Protein
Catalog #: PEH6110



Summary

Name	MS18B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O43482
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	CGGTERAIDQ ASFTTSMEWD TQVVKGSSPL GPAGLGAEPP AAGPQLPSWL QPERCAVFQC AQCHAVLADS VHLAWDLRSR LGAVVFSRVT NNVVLEAPFL VGIEGSLKGS TYNLLFCGSC GIPVGFHLYS THAALAALRG HFCLSSDKMV CYLLKTKAIV NASEMDIQNV PLSEKIAELK

image



Background

Alternative Names Protein Mis18-beta (Cancer/testis antigen 86) (CT86) (Opa-interacting protein 5) (OIP-5)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human OIP5 Protein
Catalog #: PEH6110

