


Product Name: Recombinant Human ZNF146 Protein
Catalog #: PEH6116



Summary

Name	OZF_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q15072
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	CGKVFSHKSN LTEHEHFHTR EKPFECEG KAFSQKQYVI KHQNTHTGEK LFECNECGKS FSQKENLLTH QKIHTGEKPF ECKDCGKAFI QKSNLIRHQR THTGEKPFVC KECGKTFSGK SNLTEHEKIH IGEKPFKCSE CGTAFGQKKY LIKHQNIHTG EKPYECEG

image



Background

Alternative Names Zinc finger protein OZF (Only zinc finger protein) (Zinc finger protein 146)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ZNF146 Protein
Catalog #: PEH6116

