


Summary

Name	NTM1A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BV86
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RFLREGPNKT GTSCALDCGA GIGRITKRL LPLFREVMV DITEDFLVQA KTYLGEEGKR VRNYFCCGLQ DFTPEPDSYD

image



Background

Alternative Names	N-terminal Xaa-Pro-Lys N-methyltransferase 1 (EC 2.1.1.244) (Alpha N-terminal protein methyltransferase 1A) (Methyltransferase-like protein 11A) (N-terminal RCC1 methyltransferase) (X-Pro-Lys N-terminal protein methyltransferase 1A) (NTM1A) [Cleaved into: N-terminal Xaa-Pro-Lys N-methyltransferase 1, N-terminally processed]
--------------------------	---

Note

Product Name: Recombinant Human NTMT1 Protein
Catalog #: PEH6124



For Research Use Only , Not for Diagnostic Use