


Product Name: Recombinant Human SNAPC5 Protein
Catalog #: PEH6127



Summary

Name	SNPC5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O75971
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MLSRLQELRK EEETLLRLKA ALHDQLNRLK VEELALQSMI SSRRGDEMLS SHTVPEQSHD MLVHVDNEAS INQTTLELST KSHVTEEEEE EEEEESDS

image



Background

Alternative Names	snRNA-activating protein complex subunit 5 (SNAPc subunit 5) (Small nuclear RNA-activating complex polypeptide 5) (snRNA-activating protein complex 19 kDa subunit) (SNAPc 19 kDa subunit)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use