


Product Name: Recombinant Human SPIN2A Protein
Catalog #: PEH6135



Summary

Name	SPI2A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q99865
Host	E.coli
Species	Human
Molecular Mass	30 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GQQTRAAAGR ATGSANMTKK KVSQKKQRGR PSSQPRRNIV GCRISHGWKE GDEPITQWKG TVLDQVPINP SLYLVKYDGI DCVYGLELHR DERVLSLKIL SDRVASSHIS DANLANTIIG KAVEHMFEGE HGSKDEWRGM VLAQAPIMKA WFYITYEKDP VLYMYQLLDD YKEGDLRIMP ESSESPTER EPGGVVDGLI GKHVEYTKED GSKRIGMVIH QVETKPSVYF IKFDDDFHIY

image



Background

Alternative Names Spindlin-2A (Protein DXF34) (Spindlin-like protein 2A) (SPIN-2) (SPIN-2A)

Note

Product Name: Recombinant Human SPIN2A Protein
Catalog #: PEH6135



For Research Use Only , Not for Diagnostic Use