


Product Name: Recombinant Human UBXN8 Protein
Catalog #: PEH6136



Summary

Name	UBXN8_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O00124
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SFKSPQVYLK EEEEKNEKRQ KLVRRKKQEA QGEKASRYIE NVLKPHQEMK LRKLEERFYQ MTGEAWKLSS GHKLGDEGT SQTSFETSNR EAAKSQNLPK PLTEFPSPAE QPTCKEIPDL PEEPSQTAEV VVTVALRCPS GNVLRRRFLK SYSSQVLFDW MTRIGYHISL YSLSTSFPRR PLAVEGGQSL EDIGITVDTV LILEEKEQTN

image



Background

Alternative Names UBX domain-containing protein 8 (Reproduction 8 protein) (Rep-8 protein)
(UBX domain-containing protein 6)

Note

Product Name: Recombinant Human UBXN8 Protein
Catalog #: PEH6136



For Research Use Only , Not for Diagnostic Use