


Product Name: Recombinant Human POLR2B Protein
Catalog #: PEH6137



Summary

Name	RPB2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P30876
Host	E.coli
Species	Human
Molecular Mass	13 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	YRSYKEQESK KGFDQEEVFE KPTRETCQGM RHAIYDKLDD DGLIAPGVRV SGDDVIIGKT VTLPENEDEL ESTNRRYTKR DCSTFLRTSE TGIVDQVMVT

image



Background

Alternative Names	DNA-directed RNA polymerase II subunit RPB2 (EC 2.7.7.6) (DNA-directed RNA polymerase II 140 kDa polypeptide) (DNA-directed RNA polymerase II subunit B) (RNA polymerase II subunit 2) (RNA polymerase II subunit B2)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use