

Product Name: Recombinant Human RHBDL3 Protein
Catalog #: PEH6140



Summary

Name	RHBL3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P58872
Host	E.coli
Species	Human
Molecular Mass	19 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VAACAEAERI EELEPEAEER LPAAPEDHWK VLFDQFDPGN TGYISTGKFR SLLESHSSKL DPHKREVLLA LADSHADGQI GYQDFVSLMS NKRSNSFRQA ILQGNRRLLS KALLEEKGLS LSQRLIRHVA YETLPREIDR

image

SDS

Background

Alternative Names	Rhomboid-related protein 3 (EC 3.4.21.105) (Ventrhoid transmembrane protein)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use