


**Product Name: Recombinant Human PCDHB8 Protein**  
**Catalog #: PEH6151**



## Summary

|                              |   |
|------------------------------|---|
| <b>Name</b>                  | PCDB8_HUMAN   |
| <b>Purity</b>                | Greater than 85% as determined by reducing SDS-PAGE   |
| <b>Endotoxin level</b>       | Please contact with the lab for this information  |
| <b>Tag</b>                   | C-His   |
| <b>Accession #</b>           | Q9UN66  |
| <b>Host</b>                  | E.coli  |
| <b>Species</b>               | Human   |
| <b>Molecular Mass</b>        | 13 kDa  |
| <b>Buffer</b>                | Phosphate buffered saline   |
| <b>Shipping</b>              | The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.   |
| <b>Stability&amp;Storage</b> | Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| <b>Sequence</b>              | VKVESPPG TAFPLKNAED LDIGQNNIEN YIISPNSYFR VLTRKRS DGR<br>KYPELVLDKA LDREEAELR LTLTALDGGG PPRSGTAQVY IEVVDVNDNA  |

## image



## Background

**Alternative Names** Protocadherin beta-8 (PCDH-beta-8) (Protocadherin-31)

## Note

For Research Use Only , Not for Diagnostic Use