


Product Name: Recombinant Human MRC2 Protein
Catalog #: PEH6153



Summary

Name	MRC2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q9UBG0
Host	E.coli
Species	Human
Molecular Mass	43 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DTQLDWICKI PRGTDVREPD DSPQGRREWL RFQEA EYKFF EHHSTWAQAAQ RICTWFQAEL TSVHSQAELD FLSHNLQKFS RAQEQHWWIG LHTSESDGRF RWTGDSIINF ISWAPGKPRP VGKDKKCVYM TASREDWGDQ RCLTALPYIC KRSNVTKETQ PPDLPPTALG GCPSDWIQFL NKCFQVQGQE PQSRVKWSEA QFSCEQQEAQ LVTITNPLEQ AFITASLPNV TFDLWIGLHA SQRDFQWVEQ

image



Background

Alternative Names	C-type mannose receptor 2 (C-type lectin domain family 13 member E) (Endocytic receptor 180) (Macrophage mannose receptor 2) (Urokinase-type plasminogen activator receptor-associated protein) (UPAR-associated protein) (Urokinase receptor-associated protein) (CD antigen CD280)
--------------------------	--

Product Name: Recombinant Human MRC2 Protein
Catalog #: PEH6153



Note

For Research Use Only , Not for Diagnostic Use