


Product Name: Recombinant Human SPAG11B Protein
Catalog #: PEH6158



Summary

Name	SG11B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q08648
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RHVNHSATEA LGELRERAPG QGTNGFQLLR HAVKRDLLPP RTPPYQVHIS HREARGPSFR ICVDFLGPRW ARGCTGN

image



Background

Alternative Names	Sperm-associated antigen 11B (Human epididymis-specific protein 2) (He2) (Protein EP2) (Sperm antigen HE2)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use