


Product Name: Recombinant Human STOML1 Protein
Catalog #: PEH6161



Summary

Name	STML1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UBI4
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DLNTATRMTA QNAMTKALLK RPLREIQMEK LKISDQLLLE INDVTRAWGL EVDRLVLA VE AVLQPPQDSP AGPNLDSTLQ QLALHFLGGS MNSMAGGAPS PGPADTVEMV SEVEPPAPQV

image



Background

Alternative Names Stomatin-like protein 1 (SLP-1) (EPB72-like protein 1) (Protein unc-24 homolog) (Stomatin-related protein) (STORP)

Note

For Research Use Only , Not for Diagnostic Use