

**Product Name: Recombinant Human NR5A1 Protein**  
**Catalog #: PEH6166**



## Summary

<b>Name</b>	STF1_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q13285
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	44 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	EAVRADRMRG GRNKFGPMYK RDRALKQKK AQIRANGFKL ETGPPMGVPP PPPPAPDYVL PPSLHGPEPK GLAAGPPAGP LGDFGAPALP MAVPGAHGPL AGYLYPAFPG RAIKSEYPEP YASPPQGLP YGYPEPFSGG PNVPELILQL LQLEPDEDQV RARILGCLQE PTKSRPDQPA AFGLLCRMAD QTFISIVDWA RRCMVFKELE VADQMTLLQN CWSSELLVFDH IYRQVQHGE GSILLVTGQE VELTTVATQA GSLLHSLVLR AQELVLQLLA LQLDRQEFVC LKFIILFSLD LKFLNNHILV KDAQEKANAA LLDYTLCHYP HCGDKFQQL LCLVEVRALS MQAKELYHK HLGNEPRNN LLIEMLQAKQ

## image



## Background

**Alternative Names** Steroidogenic factor 1 (SF-1) (STF-1) (hSF-1) (Adrenal 4-binding protein)

**Product Name: Recombinant Human NR5A1 Protein**  
**Catalog #: PEH6166**



---

(Fushi tarazu factor homolog 1) (Nuclear receptor subfamily 5 group A member 1) (Steroid hormone receptor Ad4BP)

**Note**

For Research Use Only , Not for Diagnostic Use