


**Product Name: Recombinant Human ST14 Protein**  
**Catalog #: PEH6167**



## Summary

<b>Name</b>	ST14_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	N-SUMO & C-His
<b>Accession #</b>	Q9Y5Y6
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	40 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	DNSCSFLGHA RGVELMRFTT PGFPDSPYPA HARCQWALRG DADSVLSLTF RSFDLASCDE RGSDLVTVYN TLSPMEPHAL VQLCGTYPPS YNLTFHSSQN VLLITLITNT ERRHPGFEAT FFQLPRMSSC GGRLRKAQGT FNSPYYPGHY PPNIDCTWNI EVPNNQHVKV RFKFFYLLEP GVPAGTCPKD YVEINGEKYC GERSQFVVTS NSNKITVRFH SDQSYTDTGF

## image



## Background

<b>Alternative Names</b>	Suppressor of tumorigenicity 14 protein (EC 3.4.21.109) (Matriptase) (Membrane-type serine protease 1) (MT-SP1) (Prostamin) (Serine protease 14) (Serine protease TADG-15) (Tumor-associated differentially-expressed gene 15 protein)
--------------------------	--

**Product Name: Recombinant Human ST14 Protein**  
**Catalog #: PEH6167**

---



**Note**

For Research Use Only , Not for Diagnostic Use