

Product Name: Recombinant Human SLC7A14 Protein
Catalog #: PEH6170



Summary

Name	S7A14_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8TBB6
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	YQPESDIDGF VKFLSEEHTK KKEGILADCE KEACSPVSEG DEFSGPATNT CGAKNLPSLG DNEMLIGKSD KSTYNVNHPN YGTVDMTTGI EADESENIYL IKLKKLIGPH YYTMRIRLGL PGKMDRPTAA

image

SDS

Background

Alternative Names Probable cationic amino acid transporter (Solute carrier family 7 member 14)

Note

For Research Use Only , Not for Diagnostic Use