


Product Name: Recombinant Human SLC05A1 Protein
Catalog #: PEH6190



Summary

Name	SO5A1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H2Y9
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RQREFPLSTV SERVGHDPNA RTRSCPAFST QGEFHEETGL QKGIQCAAQT YPGPFPEAIS SSADPGLEES

image



Background

Alternative Names	Solute carrier organic anion transporter family member 5A1 (Organic anion transporter polypeptide-related protein 4) (OATP-RP4) (OATPRP4) (Solute carrier family 21 member 15)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use