

Product Name: Recombinant Human AGPAT4 Protein
Catalog #: PEH6195



Summary

Name	PLCD_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NRZ5
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DRKTVATSLQ HLRDYPEKYF FLIHCEGTRF TEKKHEISMQ VARAKGLPRL KHLLPRTKG FAITVRSLRN VVSAVDCTL NFRNNENPTL LGVLNGKKYH ADLYVRRIP L EDIPEDDDEC SAWLHKLYQE

image

SDS

Background

Alternative Names 1-acyl-sn-glycerol-3-phosphate acyltransferase delta (EC 2.3.1.51) (1-acylglycerol-3-phosphate O-acyltransferase 4) (1-AGP acyltransferase 4) (1-AGPAT 4) (Lysophosphatidic acid acyltransferase delta) (LPAAT-delta)

Note

For Research Use Only , Not for Diagnostic Use