


Product Name: Recombinant Human SSC4D Protein
Catalog #: PEH6198



Summary

Name	SRB4D_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q8WTU2
Host	E.coli
Species	Human
Molecular Mass	27 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VGGPSRCRGR LEVMHGGSWG SVCDDWVDV DANVVCRLG CGLALPVRP LAFGQGRGPI LLDNVECRGQ EAALSECGR GWGVHNCFHY EDVAVLCDEF LPTQPPTKMLTSRAPPTTL

image



Background

Alternative Names	Scavenger receptor cysteine-rich domain-containing group B protein (Four scavenger receptor cysteine-rich domains-containing protein) (S4D-SRCRB)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use