



Product Name: Recombinant Human RFNG Protein
Catalog #: PEH6204



Summary

Name	RFNG_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q9Y644
Host	E.coli
Species	Human
Molecular Mass	30 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PAPRAPPSRP AAPSLRPDDV FIAVKTTTRKN HGPRLRLLLR TWISRARQQT FIFTDGDDPE LELQGGDRV I NTNCSAVRTR QALCCKMSVE YDKFIESGRK WFCHVDDDNY VNARSLHLHLL SSFSPSQDVY LGRPSLDHPI EATERVQGG R

image



Background

Alternative Names Beta-1,3-N-acetylglucosaminyltransferase radical fringe (EC 2.4.1.222) (O-fucosylpeptide 3-beta-N-acetylglucosaminyltransferase)

Note

For Research Use Only , Not for Diagnostic Use