


Product Name: Recombinant Human VSIG2 Protein
Catalog #: PEH6210



Summary

Name	VSIG2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q96IQ7
Host	E.coli
Species	Human
Molecular Mass	36 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EPLSTPLGKT AELTCTYSTS VGDSFALEWS FVQPGKPISE SHPILYFTNG HLYPTGSKSK RVSLQNPPT VGVATLKLTD VHPSDTGTYL CQVNNPPDFY TNGLGLINLT VLVPPSNPLC SQSGQTSVGG STALRCSSSE GAPKPVYNWV RLGTFPTPSP GSMVQDEVSG QLILTNSLT SSGTYRCVAT NQMGSASCEL TLSVTEPSQG

image



Background

Alternative Names V-set and immunoglobulin domain-containing protein 2 (Cortical thymocyte-like protein) (CT-like protein)

Note

Product Name: Recombinant Human VSIG2 Protein
Catalog #: PEH6210



For Research Use Only , Not for Diagnostic Use