

Product Name: Recombinant Human SNAI1 Protein
Catalog #: PEH6216



Summary

Name	SNAI1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-His
Accession #	O95863
Host	E.coli
Species	Human
Molecular Mass	19 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SDPNRKPNYS ELQDSNPEFT FQQPYDQ AHL LAAIPPPEIL NPTASLPMLI WDSVLAPQAQ PIAWASLR LQ ESPRVAELTS LSEDESGKGS QPPSPSPAP SSFSTSVSS LEAEAYAFP GLGQV PKQLA QLSEAKDLQA RKA FNCKYCN

image



Background

Alternative Names Zinc finger protein SNAI1 (Protein snail homolog 1) (Protein sna)

Note

For Research Use Only , Not for Diagnostic Use